



KIMBERLY-CLARK* MIC* Percutaneous Endoscopic Gastrostomy (PEG) Feeding Tubes

- Endoscopic Placement

- Medical Grade Silicone Construction
- SECUR-LOK* External Retention Ring
- Universal and Bolus Feeding Port Connectors
- Medication Port
- Collapsible Internal Retention Bumper
- Radiopaque Stripe and Bumper
- Tubing Clamp
- Available in PUSH (OTW) and PULL Placement Methods
- Ethylene Oxide (ETO) Sterilized


PUSH PEG Kit Includes:

- 1- Traction Removable PUSH PEG Tube
- 1- 5 Fr Retrieval Snare (1.8mm X 260cm)
- 1- Luer Lock Feeding Adapter with cap
- 1- Bolus Luer Lock Feeding Adapter with cap*
- 1- SECUR-Bolus Feeding-LOK* Retention Adapter* Ring
- 1- 19ga SECURX-1 LOK*.5in. Filter Retention Needle Ring
- 1- 2519ga X 15in.. Needle. Filter Needle
- 1- Fenestrated 25ga X 1in. Needle Drape with Adhesive Tabs
- 1- 12ml Fenestrated Syringe Drape with Adhesive Tabs
- 1- Scalpel 12ml Syringe with #11 Blade
- 1- Staincalpesslwith Steel#11 IntroBladucer Needle
- 1- Guidewire Stainless Steel Introducer Needle
- 1- Tubing Guidewire Clamp
- 1- Water Tubing Soluble Clamp Lubricant
- 41- 4in Water. X 4in Soluble. Gauze Lubricant Sponges
- 24- 24in. X 24in. Slit Gauze Gauze Sponges Sponges
- 12- Suture 2in. X 2in Strand. Slit Gauze Sponges
- 1- Directions Suture Strand for Use
- 1- Directions for Use

PULL PEG Kit Includes:

- 1- Traction Removable PULL PEG Tube*
- 1- 5 Fr Retrieval Snare (1.8mm X 260cm)
- (Not available in serie #0645)
- 1- Luer Lock Feeding Adapter with cap*
- 1- Bolus Feeding Adapter*
- 1- SECUR-LOK* Retention Ring*
- 1- 19ga X 1.5in. Filter Needle
- 1- 25ga X 1in. Needle
- 1- Fenestrated Drape with Adhesive Tabs
- 1- 12ml Syringe
- 1- Scalpel with #11 Blade*
- 1- Introducer Needle*
- 1- Looped Placement Wire*
- 1- Tubing Clamp*
- 1- Water Soluble Lubricant
- 4- 4in. X 4in. Gauze Sponges
- 2- 2in. X 2in. Slit Gauze Sponges
- 1- Suture Strand
- 1- Directions for Use*

Serie # 0644 contains only the products indicated with *



Item #	Size		
CE marked	(Outer Diameter)	Description	Packaging
0630-14	14 Fr	PUSH (Over-The-Wire) Method	1 Each
0630-20	20 Fr	PUSH (Over-The-Wire) Method	1 Each
0630-24	24 Fr	PUSH (Over-The-Wire) Method	1 Each

Item #	Size		
CE marked	(Outer Diameter)	Description	Packaging
0640-14	14 Fr	PULL Method	1 Each
0640-20	20 Fr	PULL Method	1 Each
0640-24	24 Fr	PULL Method	1 Each
0644-14	14 Fr	PULL Method (limited Kit)	1 Each
0644-20	20 Fr	PULL Method (limited Kit)	1 Each
0644-24	24 Fr	PULL Method (limited Kit)	1 Each
0645-14	14 Fr	PULL Method	1 Each
0645-20	20 Fr	PULL Method	1 Each
0645-24	24 Fr	PULL Method	1 Each